# **Product Specification**

**ProSpore**®

# **Product Description**

ProSpore® is a sealed, glass ampoule containing a prepared media of Soybean-Casein Digest Broth and Bromocresol Purple pH indicator which is inoculated with a sterilization-resistant spore population. When the ampoule is incubated at the appropriate temperature, the organism will grow. The physical characteristic showing growth is a color change from purple to or toward yellow and/or turbidity.

#### **Technical Parameters**

Organism	Geobacillus
e. gamen	stearothermophilus
	7953 <sup>1</sup>
Spore Count	≥1.0 x 10 <sup>4</sup>
	≥1.0 x 10 <sup>5</sup>
	≥1.0 x 10 <sup>6</sup>
D-Value	≥1.5 minutes
Z-Value	≥6 °C
Incubation	Time: 48 hrs.
	Temperature: 55-60°C
Storage	Refrigerate: 2-8°C
Shelf-Life	18 months from date of
	manufacture
Nominal	Diameter: 14.75mm
Dimensions	Height: 55mm

## **Indications for Use**

ProSpore Ampoule is a self-contained biological indicator for steam sterilization ideal for validating/monitoring the sterilization of liquid loads.

## Packaging

ProSpore is contained in a glass ampoule. These glass ampoules are packaged in boxes of 50 or 10. The boxes are made of thick, dual-surfaced chipboard. These boxes have partitions that divide them into equal divisions.

#### Labeling

- Each ampoule is labeled with product type, organism, and lot number.
- Each box contains a self-adhesive label on the outside to further identify the product.

### **Reorder Information**

- PS-4-10 (10<sup>4</sup>, 10/box)
- PS-4-50 (10<sup>4</sup>, 50/box)
- PS-5-10 (10<sup>5</sup>, 10/box)
- PS-5-50 (10<sup>5</sup>, 50/box)
- PS-6-50 (10<sup>6</sup>, 50/box)



<sup>11</sup> or equivalent recognized culture collection strain