

Product Specification

Tryptic Soy Broth with Bromocresol Purple

Product Description

TSB-BP (Mesa modified Tryptic Soy Broth with Bromocresol Purple) media is specially formulated soybean casein digest culture medium containing Bromocresol Purple as a pH indicator. The TSB-BP culture medium is specifically formulated for use with Mesa's Spore Strip brand of products and supplied in individual tubes for easy aseptic transfer and culturing. Test results are obtained during incubation by observing the medium which changes color from purple to or toward yellow and/or turbid if any spores have survived the process.

Technical Parameters

Incubation	24 hours for Steam 72 hours for Chemiclave 7 days for EO / Dry Heat
Storage	2-30°C Protect from light.
Shelf-Life	12 months from date of manufacture
Nominal Dimensions	Media Tube: 16 x 100mm Volume: 6.5mL

Indications for Use

TSB-BP medium is specially formulated for use with appropriate Mesa Labs spore tests.

Packaging

TSB-BP is packaged as 100 tubes of culture media in a chipboard box.

Labeling

- Each media tube is labeled with lot number and expiration date.
- Each box has a self-adhesive label on the outside to further identify the product.

Reorder Information

- TSB-BP16 (100 sets/box)

